

Installation Requirements:

- Matlab (v 2011 or higher)
- Java Runtime Environment (v 1.7 or higher)

1. Download the "Cropping and Visualization Tools" folder
 - a. Open "CustomCrop_1.m"
 - i. Edit input file name "fname1"
 - b. Run "CustomCrop_1.m"
 - i. It should generate an image of the input file
 - ii. Draw a line around the area of interest
 - Choose a single cell to draw around
 - This generates the mask that is used in CustomCrop_2.m
 - c. Open "CustomCrop_2.m"
 - i. Edit output file name "fname2"
 - d. Run "CustomCrop_2.m"
 - i. This generates output images of the area of interest a blurred edge
 - The blurring can be made sharper or smoother by editing the "blurmask" variable
 - e. Open "AntibodyShadowCrop.m" (to generate the same crop on separate antibody images)
 - i. Edit input file name "fname3"
 - ii. Edit output file name "fname4"
 - f. Run "AntibodyShadowCrop.m"
2. Download the "Cell Image Analysis" folder (separate user guide included in the folder for more details)
 - a. Create a folder inside the "data" folder with appropriate experiment name (eg. "test")
 - b. Specify the path of Matlab in the "config" file
 - c. Can create multiple such folders and run together
 - d. Create a folder named "actin" inside the data>"test" folder
 - e. Put all the actin images in the "actin" folder
 - i. Use files generated as output files from Step 1
 - ii. Preferable format: .tif
 - f. Open the Jar file "FilamentsAnalysis_v0.9_debug"
 - i. Click on the "Add" button
 - This will open a separate GUI window
 - Write the name of the condition/experiment
 - Click the button next to the "Filaments Image"
 - Choose the proper "actin" folder
 - ii. Multiple conditions can be added one after another

- iii. Choose the parameters in the GUI (example screenshot provided)
 - g. Click "Run All"
 - i. Progress of the analysis can be monitored
 - h. Output files will be stored on a folder named "res_filaments" inside the "actin" folder
 - i. The statistical information regarding single filaments is stored in the file named "dataM-[filename].tif.mat"
- 3. Download the "OOP Tools" folder
- 4. Preparing Angle List
 - a. Open "get_angles_for_OOP.m"
 - i. Edit input file name "fname1"
 - The input name should be same as the statistical output file from the Cell Image Analysis package ("dataM-[filename].tif.mat")
 - b. Edit output file name "fname2"
 - c. Run "get_angles_for_OOP.m"
 - d. If a batch of output files need to be processed together, use the "get_angles_for_OOP_batch.m"
 - e. Save all the output files in a specific folder (eg. "angles")
- 5. Calculate Orientational Order Parameter
 - a. Open "OOP_maincode.m"
 - b. Run "OOP_maincode.m"
 - c. Choose the folder in which output files from Step 5 were saved ("angles")
 - d. Enter the number of different conditions that will be analyzed
 - e. Enter the name of the first condition
 - f. Choose all the files (in the "angles" folder) that belong to the same condition
 - g. Repeat for all the conditions
 - h. Choose a name for the output/summary file
 - i. It will create a file "[summaryname].mat"
 - j. Open "[summaryname].mat"
 - k. Average OOP values for individual conditions are stored in "OOP_cond" variable
 - l. Standard deviation values for individual conditions are stored in "OOP_E" variable
 - m. Each column in the OOP_VEC corresponds to the OOP values of individual cells in that condition
 - i. The third column will contain OOP values of individual cells in condition #3
 - n. Number of cells analyzed in each condition is stored in the variable "num_cs"